



Product specification

Tube, 7 ml, (LxØ): 50 x 16 mm, PP, with print



Product information

Product number	62.536.009
----------------	------------

Product descriptions

Name	Tube, 7 ml, (LxØ): 50 x 16 mm, PP, with print
Product number	62.536.009
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	16 mm
Length of product	50 mm
Height of product	50 mm

Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	white

Graduation	yes
Closure type	Push cap
Autoclavable (max. °C)	121 °C

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	72000
Volume of case	0.0492 cbm
Weight of case	5.1 kg
EAN inner box	4038917110037
EAN case	4038917033107

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT