



Product specification

Tube, 12 ml, (LxØ): 110 x 17 mm, PP



Product information

Product number	57.527.300
----------------	------------

Product descriptions

Name	Tube, 12 ml, (LxØ): 110 x 17 mm, PP
Product number	57.527.300
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PP, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	17 mm
Length of product	110 mm
Height of product	110 mm

Product characteristics

Base shape	Conical base
Graduation	yes
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	3000 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0492 cbm
Weight of case	5.48 kg
EAN inner box	4038917367431
EAN case	4038917367424

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT