



## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**



### Product information

Product number	55.526.535
----------------	------------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PP
Product number	55.526.535
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, red, push cap, 1,000 piece(s)/bag

### Material & colours

Color of product	red
------------------	-----

### Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	64000
Volume of case	0.0813 cbm
Weight of case	10.3 kg
EAN inner box	4038917109857
EAN case	4038917032353

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT