



## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**



### Product information

Product number	55.526.526
----------------	------------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PP
Product number	55.526.526
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, purple, push cap, 500 piece(s)/StackPack

### Material & colours

Product material	Polypropylene (PP)
Color of product	purple

### Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	4000
Piece(s) / pallet	120000
Volume of case	0.046 cbm
Weight of case	10.252 kg
EAN inner box	4038917109802
EAN case	4038917032223

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT