



Product specification

Tube, 5 ml, (LxØ): 75 x 12 mm, PP, with print



Product information

Product number	62.526.003
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PP, with print
Product number	62.526.003
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white

Graduation	yes
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	1000000
Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	64000
Volume of case	0.0813 cbm
Weight of case	10.36 kg
EAN inner box	4038917109291
EAN case	4038917057776

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT