




## Product specification

**Carrier tube, 5 ml, (LxØ): 75 x 13 mm, PP**



### Product information

Product number	55.525.302
----------------	------------

### Product descriptions

Name	Carrier tube, 5 ml, (LxØ): 75 x 13 mm, PP
Product number	55.525.302
Product description	Carrier tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, brown, push cap, 250 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	brown

### Size

Diameter	13 mm
Length of product	75 mm
Height of product	75 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C

Centrifugation max (RCF)	5000 x g
--------------------------	----------

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	250
Piece(s) / inner box	250
Piece(s) / case	2000
Piece(s) / pallet	60000
Volume of case	0.046 cbm
Weight of case	5.72 kg
EAN inner box	4038917367257
EAN case	4038917367240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT