



Product specification

Tube, 5 ml, (LxØ): 75 x 13 mm, PP



Product information

Product number	55.525.061
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 13 mm, PP
------	-----------------------------------

Product number	55.525.061
----------------	------------

Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, with graduations, 100 piece(s)/bag
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	13 mm
----------	-------

Length of product	75 mm
-------------------	-------

Height of product	75 mm
-------------------	-------

Product characteristics

Base shape	Round base
------------	------------

Graduation	yes
------------	-----

Closure type	Push cap
--------------	----------

Autoclavable (max. °C)	121 °C
------------------------	--------

Centrifugation max (RCF)	5000 x g
--------------------------	----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	4000
Piece(s) / pallet	48000
Volume of case	0.1237 cbm
Weight of case	5.852 kg
EAN inner box	4038917276825
EAN case	4038917276818

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT