



Product specification

Tube, 5 ml, (LxØ): 75 x 13 mm, PP



Product information

Product number	55.525.060
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 13 mm, PP
Product number	55.525.060
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, without cap, with graduations, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	13 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	yes
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	50000
Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	64000
Volume of case	0.0813 cbm
Weight of case	11.6 kg
EAN inner box	4038917260008
EAN case	4038917556460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT