



## Product specification

**Tube, 6.5 ml, (LxØ): 85 x 13 mm, PP**



### Product information

Product number	55.522.300
----------------	------------

### Product descriptions

Name	Tube, 6.5 ml, (LxØ): 85 x 13 mm, PP
Product number	55.522.300
Product description	Tube, working volume: 6.5 ml, (LxØ): 85 x 13 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

### Size

Diameter	13 mm
Length of product	85 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	72000
Volume of case	0.0492 cbm
Weight of case	6.18 kg
EAN inner box	4038917367226
EAN case	4038917367219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT