



Product specification

Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP



Product information

Product number	58.505.026
----------------	------------

Product descriptions

Name	Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP
------	--------------------------------------

Product number	58.505.026
----------------	------------

Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, transparent, push cap, natural, cap enclosed, 100 piece(s)/bag
---------------------	---

Material & colours

Color of product	transparent
------------------	-------------

Color of cap	natural
--------------	---------

Size

Diameter	16.5 mm
----------	---------

Length of product	95 mm
-------------------	-------

Product characteristics

Base shape	Flat base
------------	-----------

Graduation	no
------------	----

Closure type	Push cap
--------------	----------

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched

yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	28000
Volume of case	0.0649 cbm
Weight of case	7.12 kg
EAN inner box	4038917108522
EAN case	4038917032407

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT