




## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print**



### Product information

Product number	62.476.326
----------------	------------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print
Product number	62.476.326
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: green, 1,000 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

### Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	green
Graduation	no
Closure type	Push cap

Centrifugation max (RCF)

3500 x g

## Purity & certification

Product category

no medical device | IVD

Batched

yes

## Packaging

Piece(s) / inner box

1000

Piece(s) / case

2000

Piece(s) / pallet

80000

Volume of case

0.0382 cbm

Weight of case

6.12 kg

EAN inner box

4038917368124

EAN case

4038917368117

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT