




## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print**



### Product information

Product number	62.476.178
----------------	------------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print
Product number	62.476.178
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: white, 1,000 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

### Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

### Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white
Graduation	no

Closure type	Push cap
Centrifugation max (RCF)	3500 x g

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	80000
Volume of case	0.0382 cbm
Weight of case	6.2 kg
EAN inner box	4038917107662
EAN case	4038917011846

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT