



Product specification

Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print



Product information

Product number	62.476.017
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS, with print
Product number	62.476.017
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: blue, with graduations, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	12 mm
Length of product	75 mm
Height of product	75 mm

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	blue
Graduation	yes

Closure type	Push cap
Centrifugation max (RCF)	3500 x g

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	80000
Volume of case	0.0382 cbm
Weight of case	6.2 kg
EAN inner box	4038917107358
EAN case	4038917032629

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT