



## Product specification

**Tube, 12 ml, (LxØ): 95 x 16.5 mm, PS, with print**



## Product information

Product number	62.470.300
----------------	------------

## Product descriptions

Name	Tube, 12 ml, (LxØ): 95 x 16.5 mm, PS, with print
Product number	62.470.300
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, flat base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	16.5 mm
Length of product	95 mm
Height of product	95 mm

## Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	white
Graduation	yes

Closure type	Push cap
--------------	----------

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0492 cbm
Weight of case	6.09 kg
EAN inner box	4038917368001
EAN case	4038917367998

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT