



Product specification

Tube, 8 ml, (LxØ): 100 x 13 mm, PS



Product information

Product number	55.467.300
----------------	------------

Product descriptions

Name	Tube, 8 ml, (LxØ): 100 x 13 mm, PS
Product number	55.467.300
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, 500 piece(s)/StackPack

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	13 mm
Length of product	100 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	3500 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	96000
Volume of case	0.0366 cbm
Weight of case	8.86 kg
EAN inner box	4038917367042
EAN case	4038917367035

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT