



Product specification

Tube, 8 ml, (LxØ): 100 x 13 mm, PS, with print



Product information

Product number	62.467.012
----------------	------------

Product descriptions

Name	Tube, 8 ml, (LxØ): 100 x 13 mm, PS, with print
Product number	62.467.012
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	13 mm
Length of product	100 mm
Height of product	100 mm

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white

Graduation	yes
Closure type	Push cap
Centrifugation max (RCF)	3500 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	2000
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	56000
Volume of case	0.0649 cbm
Weight of case	9.2 kg
EAN inner box	4038917106856
EAN case	4038917032599

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT