



## Product specification

**Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP**



### Product information

Product number	58.619
----------------	--------

### Product descriptions

Name	Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP
Product number	58.619
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, skirted conical base, transparent, push cap, 500 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

### Size

Diameter	16.5 mm
Length of product	95 mm

### Product characteristics

Base shape	skirted conical base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C

### Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0492 cbm
Weight of case	6.1 kg
EAN inner box	4038917082945
EAN case	4038917032421

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT