



Product specification

Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP



Product information

| | |
|----------------|--------|
| Product number | 58.619 |
|----------------|--------|

Product descriptions

| | |
|---------------------|---|
| Name | Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP |
| Product number | 58.619 |
| Product description | Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, skirted conical base, transparent, push cap, 500 piece(s)/bag |

Material & colours

| | |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 95 mm |

Product characteristics

| | |
|------------------------|----------------------|
| Base shape | skirted conical base |
| Graduation | no |
| Closure type | Push cap |
| Autoclavable (max. °C) | 121 °C |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
|------------------|-------------------------|

Batched

yes

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 500 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 36000 |
| Volume of case | 0.0492 cbm |
| Weight of case | 6.1 kg |
| EAN inner box | 4038917082945 |
| EAN case | 4038917032421 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT