



Product specification

Tube, 12 ml, (LxØ): 95 x 16.5 mm, PS



Product information

Product number	58.618
----------------	--------

Product descriptions

Name	Tube, 12 ml, (LxØ): 95 x 16.5 mm, PS
Product number	58.618
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, skirted conical base, transparent, push cap, without cap, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	16.5 mm
Length of product	95 mm

Product characteristics

Base shape	skirted conical base
Graduation	no
Closure type	Push cap

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	36000
Volume of case	0.0492 cbm
EAN inner box	4038917082938
EAN case	4038917001922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT