



Product specification

Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print



Product information

Product number	62.9924.294
----------------	-------------

Product descriptions

Name	Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print
Product number	62.9924.294
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, screw cap, yellow, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Size

Diameter	16 mm
Length of product	100 mm

Product characteristics

Base shape	Conical base
Label/Print	With print
Color of print/label	white

Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Height of inner box	10 mm
Width of inner box	380 mm
Depth of inner box	620 mm
Height of case	245 mm
Width of case	400 mm
Depth of case	400 mm
Volume of case	0.0392 cbm
Weight of case	6.2 kg
EAN inner box	4038917369954
EAN case	4038917369947

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT