



## Product specification

**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP, with print**



## Product information

Product number	62.543.030
----------------	------------

## Product descriptions

Name	Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP, with print
------	---

Product number	62.543.030
----------------	------------

Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, skirted conical base, white, screw cap, natural, cap assembled, with print, label/print: black, with graduations, 125 piece(s)/bag
---------------------	---

## Material & colours

Color of product	white
------------------	-------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	natural
--------------	---------

## Size

Diameter	25 mm
----------	-------

Length excluding cap	107 mm
----------------------	--------

Length including cap	110 mm
----------------------	--------

Length of product	107 mm
-------------------	--------

## Product characteristics

Base shape	skirted conical base
------------	----------------------

Label/Print	With print
Color of print/label	black
Graduation	yes
Light protection	yes
Closure type	Screw cap
Centrifugation max (RCF)	8800 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

## Packaging

Piece(s) / inner box	125
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	6.84 kg
EAN inner box	4038917111713
EAN case	4038917076203

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT