



Product specification

Tube, 10 ml, (LxØ): 95 x 16.8 mm, PS, with print



Product information

Product number	62.492
----------------	--------

Product descriptions

Name	Tube, 10 ml, (LxØ): 95 x 16.8 mm, PS, with print
Product number	62.492
Product description	Tube, working volume: 10 ml, (LxØ): 95 x 16.8 mm, material: PS, round base, transparent, push cap, without cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	16.8 mm
Length of product	95 mm

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white
Graduation	yes
Closure type	Push cap

Centrifugation max (RCF)

4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	36000
Height of case	310 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0491 cbm
Weight of case	6.46 kg
EAN inner box	4038917081993
EAN case	4038917001106

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT