



Product specification

Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label



Product information

Product number	63.9922.254
----------------	-------------

Product descriptions

Name	Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label
Product number	63.9922.254
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: white/yellow, Patient information, with graduations, sterile, 500 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Size

Diameter	25 mm
Length of product	90 mm

Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	white/yellow
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	15000
Height of inner box	264 mm
Width of inner box	324 mm
Depth of inner box	608 mm
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	4.98 kg
EAN inner box	4038917151573
EAN case	4038917027786

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.