



## Product specification

**Screw cap tube, 25 ml, (LxØ): 54 x 27 mm, PS**



### Product information

Product number	60.9922.115
----------------	-------------

### Product descriptions

Name	Screw cap tube, 25 ml, (LxØ): 54 x 27 mm, PS
Product number	60.9922.115
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 54 x 27 mm, material: PS, flat base, transparent, screw cap, yellow, cap assembled, sterile, 50 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

### Size

Diameter	27 mm
Length of product	54 mm

### Product characteristics

Base shape	Flat base
Graduation	no

Closure type	Screw cap
--------------	-----------

## Purity & certification

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	14000
Height of case	238 mm
Width of case	373 mm
Depth of case	458 mm
Volume of case	0.0407 cbm
Weight of case	4.8 kg
EAN inner box	4038917151238
EAN case	4038917027731

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT