

Product specification

Screw cap tube, 25 ml, (LxØ): 54 x 27 mm, PS



Product information

Order number	60.9922.115
--------------	-------------

Product descriptions

Name	Screw cap tube, 25 ml, (LxØ): 54 x 27 mm, PS
Product number	60.9922.115
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 54 x 27 mm, material: PS, flat base, transparent, screw cap, HD-PE, cap: yellow, cap assembled, sterile, 50 piece(s)/bag

Product characteristics

Base shape	Flat base
Graduation	no
Type of cap	Screw cap

Size

Diameter	27 mm
Length of product	54 mm

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Purity & certification

Satisfies the requirement	ADR
Product category	no medical device IVD
Purity standard	sterile

Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	14000
Volume of case	0.0454 cbm
Weight of case	4.8 kg
EAN inner box	4038917151238
EAN case	4038917027731

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

