



Product specification

Screw cap tube, 120 ml, (LxØ): 114 x 44 mm, PP



Product information

Product number	60.597.001
----------------	------------

Product descriptions

Name	Screw cap tube, 120 ml, (LxØ): 114 x 44 mm, PP
Product number	60.597.001
Product description	Screw cap tube, working volume: 120 ml, (LxØ): 114 x 44 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, sterile, 25 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	44 mm
Length of product	114 mm

Product characteristics

Base shape	Flat base
Graduation	no

Closure type	Screw cap
--------------	-----------

Purity & certification

Satisfies the requirement	ADR
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	3600
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	2.76 kg
EAN inner box	4038917114172
EAN case	4038917014212

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT