



## Product specification

**Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP**



### Product information

Product number	60.541.109
----------------	------------

### Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP
Product number	60.541.109
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap enclosed, sterile, 500 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Tube material	Polypropylene (PP)
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

### Size

Diameter	16.5 mm
Length of product	101 mm

### Product characteristics

Base shape	Round base with skirted base
Graduation	no
Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	28000
Volume of case	0.065 cbm
Weight of case	7.2 kg
EAN inner box	4038917110884
EAN case	4038917028806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT