



## Product specification

**Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP**



### Product information

Product number	62.732.519
----------------	------------

### Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP
Product number	62.732.519
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 76 x 20 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, label/print: white, with graduations, 100 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Size

Diameter	20 mm
Length of product	76 mm

### Product characteristics

Base shape	skirted conical base
Color of print/label	white

Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	18000
Height of case	265 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.042 cbm
Weight of case	3.81 kg
EAN inner box	4038917251617
EAN case	4038917251600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT