



## Product specification

**Screw cap tube, 6 ml, (LxØ): 92 x 11.5 mm, PP**



### Product information

|                |            |
|----------------|------------|
| Product number | 60.546.002 |
|----------------|------------|

### Product descriptions

|                     |  |
|---------------------|--|
| Name                | Screw cap tube, 6 ml, (LxØ): 92 x 11.5 mm, PP  |
| Product number      | 60.546.002   |
| Product description | Screw cap tube, working volume: 6 ml, (LxØ): 92 x 11.5 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag |

### Material & colours

|                  |                                   |
|------------------|-----------------------------------|
| Product material | Polypropylene (PP)                |
| Color of product | transparent                       |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap     | natural                           |

### Size

|                   |         |
|-------------------|---------|
| Diameter          | 11.5 mm |
| Length of product | 92 mm   |

### Product characteristics

|            |            |
|------------|------------|
| Base shape | Round base |
| Graduation | no         |

|                          |           |
|--------------------------|-----------|
| Closure type             | Screw cap |
| Centrifugation max (RCF) | 4000 x g  |

## Purity & certification

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA                    |
| Product category          | no medical device   IVD |
| Purity standard           | sterile                 |
| Sterilisation             | Electron irradiation    |
| Batched                   | yes                     |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 1000          |
| Piece(s) / inner box   | 500           |
| Piece(s) / case        | 1000          |
| Piece(s) / pallet      | 56000         |
| Height of case         | 125 mm        |
| Width of case          | 398 mm        |
| Depth of case          | 598 mm        |
| Volume of case         | 0.0298 cbm    |
| Weight of case         | 3.76 kg       |
| EAN inner box          | 4038917111768 |
| EAN case               | 4038917013376 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT