



## Product specification

**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP**



### Product information

Product number	60.545
----------------	--------

### Product descriptions

Name	Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP
Product number	60.545
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, conical base, transparent, screw cap, natural, cap assembled, sterile, 50 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Size

Diameter	25 mm
Length of product	107 mm

### Product characteristics

Base shape	Conical base
Graduation	no

Closure type	Screw cap
Centrifugation max (RCF)	8800 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	14000
Height of case	255 mm
Width of case	398 mm
Depth of case	588 mm
Volume of case	0.0597 cbm
Weight of case	6.31 kg
EAN inner box	4038917082389
EAN case	4038917001489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT