



## Product specification

**Tube, 23 ml, (LxØ): 75 x 23.5 mm, PS**



### Product information

Product number	58.490
----------------	--------

### Product descriptions

Name	Tube, 23 ml, (LxØ): 75 x 23.5 mm, PS
Product number	58.490
Product description	Tube, working volume: 23 ml, (LxØ): 75 x 23.5 mm, material: PS, flat base, transparent, push cap, without cap, 500 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	23.5 mm
Length of product	75 mm

### Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	1500 x g

### Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of inner box	110 mm
Width of inner box	645 mm
Depth of inner box	580 mm
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	7.15 kg
EAN inner box	4038917081979
EAN case	4038917001083

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT