



Product specification

Tube, 21 ml, (LxØ): 65 x 23.5 mm, PS



Product information

Product number	58.489
----------------	--------

Product descriptions

Name	Tube, 21 ml, (LxØ): 65 x 23.5 mm, PS
Product number	58.489
Product description	Tube, working volume: 21 ml, (LxØ): 65 x 23.5 mm, material: PS, flat base, transparent, push cap, without cap, 500 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	23.5 mm
Length of product	65 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	1500 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	1000
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	6.33 kg
EAN inner box	4038917081962
EAN case	4038917001076

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT