



Product specification

Tube, 21 ml, (LxØ): 65 x 23.5 mm, PS



Product information

Product number	58.489
----------------	--------

Product descriptions

Name	Tube, 21 ml, (LxØ): 65 x 23.5 mm, PS
Product number	58.489
Product description	Tube, working volume: 21 ml, (LxØ): 65 x 23.5 mm, material: PS, flat base, transparent, push cap, without cap, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	23.5 mm
Length of product	65 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	1500 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Packaging

Piece(s) / inner box	500
----------------------	-----

Piece(s) / case	1000
-----------------	------

Piece(s) / pallet	30000
-------------------	-------

Height of inner box	80 mm
---------------------	-------

Width of inner box	565 mm
--------------------	--------

Depth of inner box	635 mm
--------------------	--------

Height of case	340 mm
----------------	--------

Width of case	398 mm
---------------	--------

Depth of case	398 mm
---------------	--------

Volume of case	0.0539 cbm
----------------	------------

Weight of case	6.33 kg
----------------	---------

EAN inner box	4038917081962
---------------	---------------

EAN case	4038917001076
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT