



## Product specification

**Tube, 7 ml, (LxØ): 50 x 16 mm, PS**



### Product information

Product number	58.485.001
----------------	------------

### Product descriptions

Name	Tube, 7 ml, (LxØ): 50 x 16 mm, PS
Product number	58.485.001
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, natural, sterile, 100 piece(s)/bag

### Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

### Size

Diameter	16 mm
Length of product	50 mm

### Product characteristics

Base shape	Flat base
Graduation	no
Cap	Standard cap
Closure type	Push cap
Centrifugation max (RCF)	2500 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	32000
Height of case	145 mm
Width of case	420 mm
Depth of case	648 mm
Volume of case	0.0395 cbm
Weight of case	4.66 kg
EAN inner box	4038917108027
EAN case	4038917012126

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT