



## Product specification

**Tube, 4.5 ml, (LxØ): 75 x 12 mm, PP**



### Product information

Product number	57.477.500
----------------	------------

### Product descriptions

Name	Tube, 4.5 ml, (LxØ): 75 x 12 mm, PP
Product number	57.477.500
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, suitable for Pharmacia/LKB Gamma Counter and Gamma Counter Cobra from Packard, 1,000 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	12 mm
Length of product	75 mm

### Product characteristics

Base shape	Conical base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	5000 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	96000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	10.592 kg
EAN inner box	4038917107778
EAN case	4038917011945

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT