



## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PS**



## Product information

Product number	55.1579.002
----------------	-------------

## Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS
Product number	55.1579.002
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500, sterile, 50 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

## Size

Diameter	12 mm
Length of product	75 mm

## Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	56000
Height of case	125 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0298 cbm
Weight of case	3.18 kg
EAN inner box	4038917162746
EAN case	4038917162739

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT