



## Product specification

**Tube, 4 ml, (LxØ): 50 x 14 mm, PS**



### Product information

Product number	55.655
----------------	--------

### Product descriptions

Name	Tube, 4 ml, (LxØ): 50 x 14 mm, PS
Product number	55.655
Product description	Tube, working volume: 4 ml, (LxØ): 50 x 14 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

### Size

Diameter	14 mm
Length of product	50 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	5000 x g

### Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order quantity	3000
Piece(s) / inner box	1000
Piece(s) / case	3000
Piece(s) / pallet	108000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	6.78 kg
EAN inner box	4038917083157
EAN case	4038917002127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT