




Product specification

Tube, 11.5 ml, (LxØ): 100 x 15.7 mm, PP



Product information

Product number	55.539.250
----------------	------------

Product descriptions

Name	Tube, 11.5 ml, (LxØ): 100 x 15.7 mm, PP
Product number	55.539.250
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PP, round base, transparent, push cap, 250 piece(s)/StackPack

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	15.7 mm
Length of product	100 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	1000
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	48000
Height of inner box	174 mm
Width of inner box	102 mm
Depth of inner box	336 mm
Height of case	215 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0315 cbm
Weight of case	5.17 kg
EAN inner box	4038917110167
EAN case	4038917029391

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT