



Product specification

Tube, 3 ml, (LxØ): 75 x 10 mm, PP



Product information

Product number	55.528
----------------	--------

Product descriptions

Name	Tube, 3 ml, (LxØ): 75 x 10 mm, PP
Product number	55.528
Product description	Tube, working volume: 3 ml, (LxØ): 75 x 10 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	10 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	5000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	96000
Height of inner box	80 mm
Width of inner box	450 mm
Depth of inner box	420 mm
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	3.89 kg
EAN inner box	4038917082235
EAN case	4038917001328

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT