



## Product specification

**Tube, 5 ml, (LxØ): 75 x 12 mm, PP, with print**



### Product information

Product number	62.526.028
----------------	------------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PP, with print
Product number	62.526.028
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 1 piece(s)/blister

### Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

### Size

Diameter	12 mm
Length of product	75 mm

### Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white

Graduation	yes
Cap	Ventilation cap
Closure type	Push cap
Centrifugation max (RCF)	5000 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	18000
Height of inner box	170 mm
Width of inner box	110 mm
Depth of inner box	350 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	2.82 kg
EAN inner box	4038917109475
EAN case	4038917012492

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT