



## Product specification

**Tube, 30 ml, (LxØ): 95 x 24.8 mm, PP**



## Product information

Product number	55.517
----------------	--------

## Product descriptions

Name	Tube, 30 ml, (LxØ): 95 x 24.8 mm, PP
Product number	55.517
Product description	Tube, working volume: 30 ml, (LxØ): 95 x 24.8 mm, material: PP, round base, transparent, push cap, 250 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Diameter	24.8 mm
Length of product	95 mm

## Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	3500 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order quantity	500
Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	18000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	5.02 kg
EAN inner box	4038917082129
EAN case	4038917001212

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT