



Product specification

Tube, 30 ml, (LxØ): 95 x 24.8 mm, PP



Product information

Product number	55.517
----------------	--------

Product descriptions

Name	Tube, 30 ml, (LxØ): 95 x 24.8 mm, PP
------	--------------------------------------

Product number	55.517
----------------	--------

Product description	Tube, working volume: 30 ml, (LxØ): 95 x 24.8 mm, material: PP, round base, transparent, push cap, 250 piece(s)/bag
---------------------	---

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	24.8 mm
----------	---------

Length of product	95 mm
-------------------	-------

Product characteristics

Base shape	Round base
------------	------------

Graduation	no
------------	----

Closure type	Push cap
--------------	----------

Centrifugation max (RCF)	3500 x g
--------------------------	----------

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	18000
Height of inner box	90 mm
Width of inner box	515 mm
Depth of inner box	640 mm
Height of case	310 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0491 cbm
Weight of case	5.02 kg
EAN inner box	4038917082129
EAN case	4038917001212

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT