



## Product specification

**Tube, 13 ml, (LxØ): 100 x 16 mm, PP, with print**



## Product information

Product number	62.515.028
----------------	------------

## Product descriptions

Name	Tube, 13 ml, (LxØ): 100 x 16 mm, PP, with print
Product number	62.515.028
Product description	Tube, working volume: 13 ml, (LxØ): 100 x 16 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 1 piece(s)/blister

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

## Size

Diameter	16 mm
Length of product	100 mm

## Product characteristics

Base shape	Round base
Label/Print	With print

Color of print/label	white
Graduation	yes
Cap	Ventilation cap
Closure type	Push cap
Centrifugation max (RCF)	3000 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	13500
Height of inner box	171 mm
Width of inner box	151 mm
Depth of inner box	351 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	778 mm
Volume of case	0.0558 cbm
Weight of case	4.06 kg
EAN inner box	4038917233118
EAN case	4038917233088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT