



Product specification

Tube, 8.5 ml, (LxØ): 75 x 15.7 mm, PS



Product information

Product number	55.495
----------------	--------

Product descriptions

Name	Tube, 8.5 ml, (LxØ): 75 x 15.7 mm, PS
Product number	55.495
Product description	Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PS, round base, transparent, push cap, without cap, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	15.7 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Packaging

Piece(s) / inner box	1000
----------------------	------

Piece(s) / case	2000
-----------------	------

Piece(s) / pallet	60000
-------------------	-------

Height of inner box	120 mm
---------------------	--------

Width of inner box	560 mm
--------------------	--------

Depth of inner box	635 mm
--------------------	--------

Height of case	340 mm
----------------	--------

Width of case	398 mm
---------------	--------

Depth of case	398 mm
---------------	--------

Volume of case	0.0539 cbm
----------------	------------

Weight of case	9.01 kg
----------------	---------

EAN inner box	4038917082013
---------------	---------------

EAN case	4038917001120
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT