



Product specification

Tube, 3.5 ml, (LxØ): 55 x 12 mm, PS



Product information

| | |
|----------------|------------|
| Product number | 55.484.005 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Tube, 3.5 ml, (LxØ): 55 x 12 mm, PS |
| Product number | 55.484.005 |
| Product description | Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Color of product | transparent |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Length of product | 55 mm |

Product characteristics

| | |
|--------------------------|------------|
| Base shape | Round base |
| Graduation | no |
| Closure type | Push cap |
| Centrifugation max (RCF) | 5000 x g |

Purity & certification

| | |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 4000 |
| Piece(s) / inner box | 1000 |
| Piece(s) / case | 4000 |
| Piece(s) / pallet | 120000 |
| Height of case | 340 mm |
| Width of case | 398 mm |
| Depth of case | 398 mm |
| Volume of case | 0.0539 cbm |
| Weight of case | 8.48 kg |
| EAN inner box | 4038917107853 |
| EAN case | 4038917012027 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT