




## Product specification

**Tube, 8 ml, (LxØ): 84 x 14.5 mm, PS**



### Product information

Product number	55.473
----------------	--------

### Product descriptions

Name	Tube, 8 ml, (LxØ): 84 x 14.5 mm, PS
Product number	55.473
Product description	Tube, working volume: 8 ml, (LxØ): 84 x 14.5 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	14.5 mm
Length of product	84 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4500 x g

### Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	60000
Height of inner box	60 mm
Width of inner box	595 mm
Depth of inner box	655 mm
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	9.1 kg
EAN inner box	4038917081818
EAN case	4038917000925

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT