




Product specification

Tube, 16.5 ml, (LxØ): 130 x 17 mm, PS



Product information

Product number	55.460
----------------	--------

Product descriptions

Name	Tube, 16.5 ml, (LxØ): 130 x 17 mm, PS
Product number	55.460
Product description	Tube, working volume: 16.5 ml, (LxØ): 130 x 17 mm, material: PS, round base, transparent, push cap, 250 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	17 mm
Length of product	130 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	3500 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of inner box	50 mm
Width of inner box	595 mm
Depth of inner box	663 mm
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	9.35 kg
EAN inner box	4038917081689
EAN case	4038917000796

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT