



## Product specification

**Salivette® Cortisol, with synthetic swab, cap: blue, with paper label**



## Product information

Product number	51.1534.500
----------------	-------------

## Product descriptions

Name	Salivette® Cortisol, with synthetic swab, cap: blue, with paper label
------	-----------------------------------------------------------------------

Product number	51.1534.500
----------------	-------------

Product description	Salivette® Cortisol, with synthetic swab, cap: blue, with paper label, label/print: white/blue, for cortisol determination from saliva, 100 piece(s)/bag
---------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------

## Material & colours

Product material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	blue

## Size

Diameter	16.8 mm
Length excluding cap	84 mm
Length including cap	97 mm

Design	with synthetic swab
Label/Print	With paper label
Color of print/label	white/blue
Closure type	Push cap
Centrifugation recommended	1000 x g
Centrifugation time	2 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	500
Type of inner packaging	StackPack
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	28000
Height of inner box	119 mm
Width of inner box	101 mm
Depth of inner box	226 mm
Height of case	140 mm
Width of case	263 mm
Depth of case	533 mm
Volume of case	0.0196 cbm
Weight of case	4.09 kg
EAN inner box	4038917159364
EAN case	4038917159357

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT