




Product specification

ESR pipette, with print



Product information

Product number	86.1996.062
----------------	-------------

Product descriptions

Name	ESR pipette, with print
Product number	86.1996.062
Product description	ESR pipette, with printed graduation and assembled piston insert, with print

Product characteristics

Label/Print	With print
-------------	------------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	42000
Height of inner box	76 mm
Width of inner box	164 mm
Depth of inner box	233 mm

Height of case	250 mm
Width of case	338 mm
Depth of case	388 mm
Volume of case	0.0328 cbm
Weight of case	8.61 kg
EAN inner box	4038917135047
EAN case	4038917025966

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT