




# Product specification

## Capillary holder



### Product information

Product number	95.1048
----------------	---------

### Product descriptions

Name	Capillary holder
Product number	95.1048
Product description	Capillary holder, capillary blood collection, 5 piece(s)/bag

### Material & colours

Color of product	red
------------------	-----

### Product characteristics

Type of collection	capillary
--------------------	-----------

### Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

### Packaging

Piece(s) / case	1000
Piece(s) / pallet	108000
Height of case	195 mm
Width of case	223 mm

Depth of case	308 mm
Volume of case	0.0134 cbm
Weight of case	2.31 kg
EAN inner box	4038917085120
EAN case	4038917003810

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT